Skin staple remover

ZILLION[®]
BLACK


Removedor de agrafos **REF 4891-521**

FR Ôte-agrafes cutanées
DE Hautklammerentferner
NL Huidnietjesverwijderaar
IT Rimuovi graffe
ES Quitagrapas
SV Agraffborttagare
BG Инструмент за отстраняване на хирургически скоби
HR Skidač kožnih spajalica
CS Odstraňovač kožních svorek

ZILLION[®] 1) **UKRP** 3) **LDPE** 4)

Bastos Viegas, s.a. IFU551999117, Rev.0 (28-11-2024)
Avenida da Fábrica, 298. 4560-164 Guilhufe, Penafiel - Portugal
Tel.: 255 729 500 - geral@bastosviegas.com - www.bastosviegas.com

EN Fig1-Exact position of the staple before removal
PT Fig1-Posição exacta do agrafó antes da remoção
FR Fig1-Position exacte de l'agrafe avant ablation
DE Abb.1:Genaue Position der Klammer vor dem Entfernen
NL Afb.1-Exacte stand van het nietje vóór verwijdering
IT Fig1-Posizione esatta della pinza prima della rimozione
ES Fig1-Posición exacta de la grapa antes de quitarla
SV Fig1Agraffens exakta läge före borttagning
BG Фиг. 1 – Точно положение на скобата преди премахане
HR Slika 1 - Točan položaj spajalice prije njezinog uklanjanja
CS Obr. 1 - Přesná poloha svorky před jejím odstraněním

<https://www.bastosviegas.com/downloads/ifu/IFU551999117.pdf>

EN SKIN STAPLE REMOVER

INTENDED PURPOSE AND INDICATIONS

Skin staple remover is intended for removal of surgical skin staples.

INSTRUCTIONS FOR USE

- 1-Inspect the staples and make sure they are perpendicular to skin before removal.
- 2-Place the lower jaw of staple remover between staple and skin (figure 1).
- 3-Squeeze and lift the staple remover perpendicularly to the skin.
- 4-The staple will be bent open by squeezing (figure 2).
- 5-Lift the open staple straight up from the tissue and discard according to local requirements.

WARNINGS AND PRECAUTIONS

- Intended to be used exclusively by healthcare professionals.
- Visually inspect the sterile barrier system prior to use. Do not use the device if you detect any breach in packaging integrity.
- Visually inspect the device. Do not use the device if you detect any sign of deterioration or defect that may compromise the performance.
- Do not re-use. Intended to be used one single time to prevent risks from contamination, including possible infections, and/or functional damage.
- Dispose according to applicable legislation and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we must inform that serious incidents related with the use of the device, that are life-threatening or compromise health, including illness or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your Country.

LABELLING SYMBOLS

1)Trademark 2)Opening direction 3)United Kingdom authorized representative

4)Polyethylene (identification of packaging material)

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